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
METHOD FOR DETERMINING OF AMITRAZ AND ITS METABOLITES IN ANIMAL PRODUCTS AND HONEY (MVI.MN 6330-2020)

Developers: Palianskikh Alena, Belyshava Ludmila, Bulgakava Volga

Scope of application: food safety control;

Purpose: quantitative determination of amitraz and its metabolites (2,4-dimethylaniline (DMA), N-(2,4-dimethylphenyl) -N'-methylformamidine (DMPF) and 2,4-dimethylformamidine (DMF)) in animal products and honey ;

Amitraz is used to treat varroaosis in farm animals and bees. Residual amounts of amitraz and its metabolites may be present in livestock products and honey.

Maximum permissible level (MPL) for amitraz and the amount of metabolites: from 0.01 to 0.40 mg / kg, depending on the type of food product.

Chromatography conditions

Agilent 1200 Chromatograph and Agilent 6410 Mass Spectrometer Detector

Stationary phase:	chromatographic column Poroshel HPH-C18 (4.6 mm □ 100 mm), grain size 2.7 μm,
Mobile phase of the composition	A (0.1% aqueous formic acid solution), B (0.01% acetonitrile ammonia solution)
Mobile phase supply conditions	gradient elution
Mobile phase feed rate	0,5 ml/min
Column oven temperature	40 °C
Injected sample volume	10 μl

Mass spectrometric detection parameters

Analyte	"Parental" ion, m/z	"Daughter" ions, m/z	Fragment tension, B	Collision energy, B
Amitraz	294,2	163,1 122,1	25	13 34
Amitraz-D3	297,2	166,1 122,1	30	14 38
DMA	122,1	107,1 79,1	70	18 26
DMPF	163,1	122,1 107,1	85	18 26
DMF	150,1	107,1 132,1	85	23 13


Sample preparation conditions

1) 1 g honey + water + internal standard + 3 n NaOH

2) 5 g milk + internal standard + 3 n NaOH

3,4) 1 g homogenized sample (liver, kidney) + water + internal standard + acetonitrile

Types of food products: 1) honey;
2) milk;
3) kidneys, liver;
4) raw fat

The proposed conditions allow:

- **to achieve the maximum response** and satisfactory chromatographic separation of amitraz and its main metabolites DMFP, DMF, DMA by using a chromatographic column filled with a reversed-phase sorbent resistant to destruction in a wide pH range, as well as a mobile phase consisting of an aqueous solution of formic acid and acetonitrile ammonia solution in a gradient eluent supply mode;
- **to achieve the maximum extraction** of four analytes from food products and the maximum purification of the extract from interfering impurities due to the stage of liquid-liquid extraction of analytes from aqueous solutions of honey samples on cartridges filled with diatomaceous earth using dichloromethane as an extractant;
- **to achieve concentration of analytes** by distilling off the solvent and dissolving the dry residue in acetonitrile and water of a smaller volume.

The range of determined values of the mass fraction of amitraz and its metabolites is:

for milk from 0.001 to 0.012 mg / kg,

for honey, liver, kidneys, raw fat - from 0.025 to 0.250 mg / kg.

Republican Unitary Enterprise «Scientific Practical Centre Of Hygiene»

Virtual exhibition of innovative developments «Hygiene safety»

Republic of Belarus, 220012
Minsk 8, Akademicheskaya Str.
Fax: +375 17 272-33-45

Scientific organization department

+375 17 310 72 91

rspch@rspch.by
edu@rspch.by

Educational centre

+375 17 399 87 24

www.rspch.by
www.certificate.by

Detailed information about the development:

Laboratory of Food Chemistry

+375 17 379-13-80
chf@rspch.by

Information about all developments
of the Centre is available at the link:
<https://rspch.by/ru/DevelopedDocuments>